

Features of The Prometheus Group TR Series and MR Series Surface EMG Biofeedback Devices

	TR Series TR-10, TR-10C, TR-20, TR-20C	MR Series MR-10, MR-15, MR-20, MR-25
	Single Channel = 10 & Dual Channel= 20	Single Channel = 10, 15 Dual Channel = 20, 25
Microvolt Activity Display	LED Bar Numeric Display on Face of Device	LED Bar and Numeric Display
Microvolt Scale	1-800 microvolts Or 1-30 microvolts Scale is imprinted on face of Device	Auto-Scale According to muscle activity
Visual Display for Threshold (Goal)	LED Bar at the Threshold Value	LCD Bar at the Threshold Value for the MR-10 & 20 LED Bar at the Threshold Value for the MR-15 & 25
Audio Tones for Threshold (Goal)	<ul style="list-style-type: none"> ▶ Continuous ▶ Above Threshold ▶ Below Threshold 	<ul style="list-style-type: none"> ▶ Above Threshold ▶ Below Threshold
Other Feedback Options		▶ Maximum Microvolt obtained, "MAX"
Work/Rest Options	<ul style="list-style-type: none"> ▶ 5 sec.Work/10 sec.Rest ▶ 10 sec.Work/10 sec.Rest 	▶ Set your own Work/Rest cycles and # of Repititions; 3 sec.to 99 sec.
Save Sessions	▶ "C" units only	▶ Saves 32 sessions – Data available on the device (without a computer download)
Computer Interface Options	<ul style="list-style-type: none"> ▶ "C" units only; ▶ Session Data only 	<ul style="list-style-type: none"> ▶ Download saved sessions ▶ Interface with a PC and the EMG Activity is displayed on the Computer Screen

Vendor Information:

Current Technology, Inc.
 2001 Regent Avenue No.
 Suite 100
 Minneapolis, MN 55422
 Telephone:800-322-4461
 Fax:763-588-0066



TR-10 / TR-10C
Single Channel



TR-20 / TR-20C
Dual Channel



MR-10



MR-15

(Dual Channel Models: MR-20 / MR-