



Thought Technology Ltd.

Distributor: Current Technology, Inc.
 2001 Regent Avenue North, Suite 100, Minneapolis, MN 55422
 Telephone: 800-322-4461 Fax: 763-588-0066
 Email: info@currenttechnologyinc.com website: www.currenttechnologyinc.com

BioGraph Infiniti PC System Requirements:

	Minimum	Recommended
Processor	1.8 GHz, Dual-core CPU	Quad-core CPU, Intel i5
Operating systems (supported) (32 and 64 bit versions)	<ul style="list-style-type: none"> - Microsoft® Windows® 8, Pro, Enterprise - Microsoft® Windows® 7 Home Premium, Professional, Enterprise or Ultimate with SP1. - Microsoft® Windows Vista® Home Basic, Home Premium, Business, Ultimate, or Enterprise with Service Pack 2. - Microsoft® Windows® XP Prof. or Home Ed. with SP3. 	
Memory (RAM)	1 GB (Windows® XP) or 2 GB (Vista® or Windows® 7/8).	2 GB (Windows® XP) or 3 GB (Vista® or Windows® 7/8).
Hard Disk space	2.5 GB available to install and run	50 GB or more to save video and sessions data
Video / Graphics card	With 256MB memory, with 1024x768 resolution support.	At least 1 GB of dedicated memory. (Example: ATI Radeon 5770, with 1280x1024 resolution and dual monitor support.)
Monitor	At least one which supports resolution of 1024 x 768 and higher.	Two monitors which support resolution of 1280x1024 and higher.

Other necessary components:

- DVD-ROM drive.
- Adequate number of USB ports depending on the intended system usage (1-4 for Pro/Flex Infiniti encoders, 2 for webcams, 1 for TT AV Sync, mouse, keyboard, printer)
- Mouse or any other compatible pointing device.
- 32-bit compatible sound card & speakers.
- Compact Flash Reader (for use with compact flash card only).
- Up to two webcams with built-in microphone (if recording video). Frame rate 30fps, Microsoft® recommended.
- Microsoft® Office with Word® and Excel® (for report generation and printing).
Microsoft Office 'Home and Student' DOES NOT WORK with the BioGraph Infiniti reporting functions.
- High-speed Internet access (for updating software and receiving technical support).



Thought Technology Ltd.

Distributor: Current Technology, Inc.

2001 Regent Avenue North, Suite 100, Minneapolis, MN 55422

Telephone: 800-322-4461 Fax: 763-588-0066

Email: info@currenttechnologyinc.com website: www.currenttechnologyinc.com

PLEASE NOTE:

- Tele-Infiniti CF (T9600) version 3.2.1 or higher is compatible with 64-bit Windows® operating systems. Version 3.1 or lower is compatible with 32-bit Windows®.
- When using Windows® 8 or Windows® 7 Home Basic or Professional or Windows Vista® Home Basic or Business, additional codecs must be installed to view camera and DVD.
- Low speed CPUs (1.8 GB) may be sluggish when running screens with big number of instruments per screen.
- If using 2 webcams for video recording, it is strongly suggested to use recommended PC requirements.
- If using TT-AV Sync functionality, please refer to possibly additional system requirements from appropriate suite manual.

The following apply to BioGraph Infiniti 6.0 and its application suites:

- *BioGraph Infiniti is a 32-bit software program.*
- *BioGraph Infiniti software installation must run with Administrator privileges.*
- *BioGraph Infiniti software and application suites are installed for **All Users** on a per-machine basis.*
- *The BioGraph Infiniti system must be installed on the PC where the database is/will be located.*
- *BioGraph Infiniti is not designed for use in a networked environment*